Genu Direxa



Features Denefits

Narrow aluminium splints and polycentric joints

Medium stabilisation of the knee without limitation of movement

Innovative TriTech material

Very comfortable, secure fit

Two elastic straps for optional use

Optimum fit for any leg shape

Covered joint bars

Suitable for contact sports

Lateral and medial grips

Straightforward positioning

Art. no. 8356

Order no.: Article number=size | Example for ordering: 8356=L-7

Size	Lower leg circumference A		Thigh circumference B	
	(cm)	(inch)	(cm)	(inch)
XXS	29-32	11.4-12.6	36-40	14.2-15.8
XS	32-35	12.6-13.8	40-44	15.8-17.3
S	35-38	13.8-15.0	44-48	17.3-18.9
М	38-41	15.0-16.1	48-52	18.9-20.5
L	41-44	16.1-17.3	52-56	20.5-22.0
XL	44-48	17.3-18.9	56-61	22.0-24.0
XXL	48-51	18.9-20.1	61-67	24.0-26.4

Indications

Knee pain with slight to moderate ligament laxity and /or slight to moderate feeling of instability

- Posttraumatic and postoperative instability and/or feeling of instability
- Knee osteoarthritis with starting to moderate degenerative ligament instability and/or slight to moderate feeling of instability (cruciate ligament focused)
- · Rheumatoid arthritis
- Patellofemoral pain syndrome (e.g. chondropath patellae, chondromalacia patellae, retropatellar osteoarthritis)

Effects

- · Improves neuromuscular function with targeted soft tissue compression
- · Relieves and stabilises the knee joint
- Improves proprioception
- Supports sensorimotor function
- · Promotes the resorption of oedemas and haematomas
- · Can relieve pain