

BELSANA

BELSANA comfortis. The latest comfort class in compression therapy.

NEW



- Soft to the touch, highly effective
- + New comfort zones
- + More individual solutions
- = **Higher customer satisfaction**



You could never be happier recommending this compression stocking than any other.


These features will convince your customers:

- ✓ Super-soft material:
For **maximum comfort** when worn and a particularly **pleasant feeling on the skin**
- ✓ New comfort zones around the foot and knee:
Comfortable fit and **no cutting-in** when sitting
- ✓ Soft tip and comfort seam: **No pressure** and **maximum freedom for the toes!**
- ✓ Reinforced around the heel and ball of the foot for **improved durability**
- ✓ **Modern appearance** in eight fresh and fashionable colours

BELSANA comfortis

Modern compression stockings with a feel-good factor

- Standard sizes I-VI plus custom sizes
- Material: 67 % polyamide, 33 % elastane
- AF and AG with elastic cuffs
- AG/T with drawstring
- AT/H with horizontal slit



CCL 1	Standard Item Number	21 5101	-	21 5104	-	21 5141	21 5143	-	-
	Made-to measure Item Number	22 5101	22 5103	22 5104	22 5147	22 5141	22 5143	-	-
CCL 2	Standard Item Number	21 5301	-	21 5304	21 5347	21 5341	21 5343	21 5346	-
	Made-to measure Item Number	22 5301	22 5303	22 5304	22 5347	22 5341	22 5343	22 5346	21 5345

Standard version	<ul style="list-style-type: none">• closed soft tip• 1 pair• short foot (shoe size 36-40)• A-F and A-G with elastic cuffs• A-G/T with drawstring• A-T/H with horizontal slit			
Leg length	A-D	Shaft (ID):	short (34-38 cm)	long (38,1-43 cm)
	A-G and A-G/T	Crotch (IG):	short (65-74 cm)	long (74,1-84 cm)
	A-T	Shaft (IK):	short (69-78 cm)	long (78,1-88 cm)
Colours	<div><div>almond</div><div>cinnamon</div><div>brittle*</div><div>wild plum blue*</div><div>elderberry blue</div><div>truffle grey*</div><div>liquorize (black)</div><div>coconut white (undyed)</div></div>			

* up to 14 days delivery time!

Special designs for standard range and size

Open tip with comfort seam in two lengths

- Short foot length (A-Y): 15-18,4 cm
- Long foot length (A-Y): 18,5-22 cm

Long foot (shoe size 41-45)

Grip tops in AD stockings:

- Effect (6 cm) NEW
- Silicone (3 or 5 cm)

Grip tops in AF and AG stockings:

- Effect (6 cm) NEW
- Silicone (5 cm)
- Silicone „Fleur“ (5 cm)
- Delivery time of 1 week for individual items in the standard range, no returns or exchanges

Special designs for made-to-measure

Asymmetric soft tip

Closed variable soft tip (5-10 cm)

Grip tops in AF and AG stockings:

- Rubber (3 oder 5 cm)

Grip tops in AD stockings:

- Rubber (5 cm)

Panty part in AT:

- With compression
- With stoma opening
- With open crotch area
- Men's stockings with vertical slit
- Reinforced crotch area

Panty part in A-T/U:

- Open





The compression stocking developed for the needs of your customers.

NEW

for more comfort

Closed soft tip

Grip top „Effect“

Knee comfort zone

Open toe with comfort seam – available in two lengths

Knee high stocking with a practical grip top also available in the standard range

Open panty part

Extra-wide panty part

Modern compression therapy developed for:

- Advanced venous disorders
- Weak connective tissues
- Maximum coverage of individual patient requirements
- High standards in terms of appearance
- Excellent wearability – for both men and women

Individual made-to-measure solutions, greater choice in the standard range.

- ✓ For patients with toes that are sensitive to pressure, or who have hallux valgus or diabetes:
 - **Variable soft tip** with manual specification of length
 - **Asymmetric soft tip shape** to fit toes
 - **Open toe in two lengths** with a comfort seam
- ✓ For patients looking for an adhesive strip solution:
 - **Practical grip top** with four effects:
No folding over, does not wear out, secure hold, better air circulation
- ✓ For patients with particular requirements for the panty part:
 - **Extra wide or open** for pregnant women or people with a very high bodyweight
 - **Compressed** for varicose veins in the genital or groin area
 - With integrated **stoma opening** for those with a stoma
 - With an **open crotch** area for bladder weakness



Tights that keep their shape



Variable soft tip



Compressive panty part



Reinforcement around the heel and ball of the foot

3 steps to the right stocking size

- 1 Find out in which size the cB measure fits
- 2 Compare each measure of the patient with each of the other measures of the defined size
- 3 No discrepancy? Size order!
Discrepancies? Made-to-measure-order


Sizes BELSANA comfortis

Circumference (cm)

extra wide panty part	75–100 70–90	80–110 75–105	85–120 80–110	95–130 90–120	100–135 95–125	105–140 100–130
extra wide panty part	90–115 90–105	95–125 95–115	100–140 100–130	105–145 105–135	115–155 115–145	120–165 120–155
Tights A–T Plus size A–G	46–67 56,1–67 46–56	49–71 60,1–71 49–60	52–74 63,1–74 52–63	56–79 67,1–79 56–67	58–84 72,1–84 58–72	62–89 76,1–89 62–76
	39–53	41–56	45–60	46–62	47–64	49–68
	31–40	35–43	37–47	39–49	41–52	43–56
	28–35	30–38	33–41	35–44	37–47	39–51
	29–36	32–40	34–43	37–45	39–48	41–52
	23–29	25–32	27–34	30–37	33–40	37–44
	19–20	20,1–22	22,1–24,5	24,6–27	27,1–30	30,1–33
	28–32	29–34	30–35	32–36	33–38	34–39
	19–23	20–24	21–25	22–26	23–27	24–28
	I	II	III	IV	V	VI

Measuring Point

Length (cm)



open toe

short foot: 15–18,4 cm
long foot: 18,5–22 cm

closed toe

short foot: 22,5–25,9 cm
long foot: 26–29,5 cm

in compliance with the most up-to-date
German RAL requirements



BELSANA Medical Products

Laubanger 23 · D-96052 Bamberg · Tel. +49 (0)951 309303-761
Fax +49 (0)951 309303-465 · infoexport@belsana.de · www.belsana.com