


BELSANA

Export Product Range







Should you have any questions, wishes or suggestions, our BELSANA export team will be more than willing to help you. You can send us at any time orders and written messages per fax or e-mail.

We can be contacted as follows:

Phone + 49 951 30 93 03-761

Fax + 49 951 30 93 03-465

E-Mail infoexport@belsana.de


www.belsana.com

BELSANA

Contents

INFORMATION	4-13
MEDICAL COMPRESSION STOCKINGS	14-3
ACCESSORIES AND CARE	32-33
SUPPORT STOCKINGS	34-4
SPORT SOCKS AND STOCKINGS	42-44
DIABETES	4.
GENERAL TERMS AND CONDITIONS	46-47

Come and visit us in Bamberg!







Innovative products of high value made in Bamberg, Germany

With **BELSANA Medical Products** you have found a reliable and strong partner in the health care market. Since 1991 BELSANA has supported customers world-wide in providing venous patients with medical compression stockings.

Under the umbrella of INDUS Holding AG, BELSANA is one of the most prestigious specialist providers in the German pharmacy market. Right from the start, the aim was to provide a range of **medical compression stockings** focussing on the individual needs of patients with vein disorders. The range of modern compression stockings from BELSANA has meanwhile become one of the most important lines in the entire market for health aids: Fine qualities, a wide selection of colours, individual made-to-measure solutions and rapid delivery give our customers the best possible support in their endeavour to comply optimally with the requirements of compression therapy.

With **pioneering innovations**, such as compression stockings with microfibres (BELSANA microsoft), comfort zones (BELSANA comfortis), merino wool (BELSANA impuls) or contemporary support and compression stockings with moisturising aloe vera (BELSANA vital Editi-

on and BELSANA vivere), BELSANA guarantees high-tech products of the highest quality. Besides this, the basic programme has been supplemented by interesting complementary products over the years. In particular, the attractive and extensive range of over-the-counter support and **travel stockings** ensures customers are offered products of the highest quality in the area of prophylaxis as well.

And what's more, BELSANA has just the right products when it comes to sports. BELSANA sport are **sport socks** and stockings with True Medical Compression, ensuring less fatigue, faster regeneration and enhanced perfor-

BELSANA compression and support stockings are **made in** Bamberg, Germany, in line with the strictest standards of quality. The company is therefore certified according to DIN EN ISO 13485, its products bearing the CE mark. The medical compression stockings from BELSANA also correspond to the strictest international quality and test provisions of RAL-GZ 387/1 and 2. Refined quality management, coupled with the development of new top products and first-class service for our customers, means that the company can grow on a firm and lasting basis.

BELSANA Quality control

The quality of medical compression stockings contributes significantly not only to wearing comfort but also to medical efficacy. Modern technological innovation ensures the high quality level of our products as documented and certified regularly by independent testing institutes.



RAL Quality Seal

Our medical compression stockings meet the standards of the RAL-Norm, the strictest norm worldwide for the production of medical compression stockings.

The RAL Quality Seal on the packaging of our medical compression stockings guarantees that our products meet this norm. The RAL guidelines are set by the "Gütezeichengemeinschaft" who decides, in cooperation with the "German Society of Phlebology", the standards (quality and technical specifications) required for the production of medical compression stockings.

The independent Hohenstein research institute carries out tests to verify the requirements set by RAL are met and principally on the basis of:

- pressure gradation
- fibre analysis
- knitting structure
- chemical analysis

OEKO-TEX ® CONFIDENCE IN TEXTILES STANDARD 100 97.0.6907 HOHENSTEIN HTII Geprüft auf Schadstoffe. www.oeko-tex.com/standard100

The Öko-Tex Standard 100 is a testing and certification system for textile products. Öko-Tex Standard 100 guarantees that certified products observe strict limits on the use of certain materials and dyes. This certified standard applies to all BELSANA compression stockings.

((

CE compliance with European directives

The CE marking serves as a "passport" for the Single European Market and allows BELSANA products to circulate freely throughout the European market. It indicates compliance with the requirements for medical products according to the European directive 93/42 EWG and the German law on medical products (MPG).

DIN EN ISO 9001 / DIN EN ISO 13485

Our products meet the highest quality standards and are produced according to the quality management system DIN EN ISO 9001 / DIN EN ISO 13485. This certificate documents process oriented operations and the optimisation of in-house production structures. The continuous improvement of processes is guaranteed by regular external audits.

Medical compression stockings vs. support stockings

BELSANA products are divided into two categories:

Medical compression stockings and support stockings. Both product categories principally work in the same way: they help the blood in the veins to flow back from the foot to the heart by exerting pressure on the legs. The pressure is the highest in the ankle and decreases in the direction of the heart. This accelerates the return circulation to the heart and the blood circulation in the leg is significantly improved.

Medical compression stockings

- While support stockings can be purchased over the counter, medical compression stockings are generally prescribed by a physician for therapeutic purposes in the case of a venous disease.
- Medical compression stockings guarantee a precisely defined and appropriate level of pressure on the dilated veins of the legs. They should be worn on a long-term daily basis.
- On the basis of his diagnosis, the physician determines the compression class as well as the specific application. The spectrum ranges from knee high stockings to tights and special designs for men and pregnant women are
- It is essential to find the most suitable design for each individual to guarantee optimal efficacy.
- BELSANA medical compression stockings are available in standard sizes as well as made-to-measure. The correct choice is made by taking the precise measurements of the patient's legs.

Support stockings

- By contrast, support stockings are worn by healthy persons as a preventive measure against venous problems.
- They are thinner than compression stockings and the pressure on the leg is not as high as with compression stockings. Nevertheless, they actively support the blood flow in the legs.
- They are optically indistinguishable from normal stockings and are recommended especially for all activities involving prolonged sitting or standing or during longer flights or car journeys.

Anti embolism stockings (TED stockings)

Both categories should not be mistaken for antiembolism stockings (also known as TED stockings – Thrombo-Embolic Deterrent), which are designed for bed-confined patients and aim at preventing thrombo-embolic disease.



The blood is pumped to the heart. The venous valves prevent the blood from flowing back.



If the vein is weak or damaged, the venous valves do not close – the blood drains downwards.



Correctly-dosed pressure has a positive effect on the stretched vein walls and promotes blood flow to the heart.

6


Compression classes and designs

The level of pressure exerted on the legs depends on the seriousness of the patient's condition: in principle, the greater the damage to the venous system is and the softer the connective tissue has become, the stronger and thicker the compression stockings must be.

The RAL-Norm distinguishes between four compression classes (Ccl.)

Compression class	Indications						
1 = slight compression	as a preventive measure						
(18–21 mm Hg)	heavy and tired legs						
	 minor varicose veins without any significant edema susceptibility 						
	varicose veins during pregnancy						
2 = medium compression	greater discomfort						
(23–32 mm Hg)	 more pronounced varicose 	veins with an edema susceptibility					
	• post traumatic swelling						
	after the healing process of insignificant ulcerations						
	• after superficial thrombophlebitis, sclerotherapy and operations on varicose veins						
3 = firm compression	constitutional or post throm	botic venous insufficiency					
(34–46 mm Hg)	 serious edema susceptibility 	,					
	secondary varicose veins						
	white atrophy						
	after the healing process of serious, recurrent ulcerations						
4 = very firm compression	lymphatic edema						
(49 mm Hg and more)	 elephantiasis 	Instead of one stocking in compression class 4					
		it is recommended that two stockings in lower					
		compression classes can be worn one on top o					

the other, as they are easier to put on.



All indications and contra-indications can be found in our packaging inlays.



Taking measurements for round- & flatknitted compression stockings

For a perfect fit and medically proper pressure gradation, it is essential that leg measurements be taken with due care and precision. Perfectly fitting stockings prevent the development of congestion and improve the acceptance by patients of compression therapy.

- Legs must be free of clothing and free of edema (finger pressure test in the ankle area – an indentation will persist after the pressure is relieved in case of edema)
- Measurements should be taken in the morning or immediately upon the removal of a compression bandage
- The leg should be measured in a standing position
- Each leg must be measured separately
- The circumference should be measured first (mark measuring points with a suitable pen, e.g. BELSANA kohl pencil)
- Measurements should be taken from the feet
- Flatknitted compression stockings have to be measured tighter than roundknitted stockings.
 For special advice, please contact the export team

		Crotch to wais	T \
		-	Γh sTv \
Left Leg		(front) (b	pack)
Length	Circumference	1 •	
IT	cT	sTh	
		H	Right Leg
IH	сН		Circumference Length
IG	cG	G	cG IG
IF	cF		cF IF
IE	cE	F)	cE IE
ID	cD		cD ID
IC	cC		cC IC
IB1	cB1		cB1 IB1
IB	сВ	B) IK	сВ ІВ
	cY	B Inside leg	сҮ
	cA (A	left or right	cA
Left foo	t		Right foot
Foot length (A-Y)		inside leg length	Foot length (A-Y)
Whole fo	ot length	Shoe size	Whole foot length

Please note: when choosing the variable soft toe (BELSANA comfortis) please fill in the whole foot length and the foot length (A-Y) for both feet. The difference must be a minimum of 5 cm and a maximum of 10 cm!

Please see measuring points stipulated


in the table. (all measurements in cm)

Measuring point	Circumference (= c)	Length (=	= 1)		
		IT* from heel to circumference cT			
(1)	cT waist	sTv sTh	from crotch to cT (front) measured from gluteal fold to cT (back) vertically		
		IH*	from heel to circumference cH		
(H)	cH hip, circumference of the trunk at its widest point (buttocks)	IK* / inside leg	from heel to crotch measured on the leg (left or right)		
G	cG no precise stipulation, at top of thigh (2 fingers wide below the crotch)	IG	from heel to circumference cG		
F	cF no precise stipulation (approx. 15 cm above cE)	IF	from heel to circumference cF		
E	cE around the knee-cap and hollow of the knee (with knee bent slightly)	IE	from heel to circumference cE		
D	cD below the knee-cap (2 fingers wide under the hollow of the knee)	ID	from heel to circumference cD		
©	cC widest calf circumference	IC	from heel to circumference cC		
B_1	cB ₁ base of the calf (in the middle of cB and cC)	IB ₁	from heel to circumference cB ₁		
В	cB ankle (the slimmest part of the lower leg, smallest leg circumference)	IB	from heel to circumference cB		
Y	cY heel – instep region (foot should be bent upwards as far as possible)				
A	cA circumference at the base of the little toe	Foot length Whole foot lenght	from heel to circumference cA from heel to the tip of the toes		

Taking measurements for flatknitted compression foot and toe part

- To measure the toes your patient should be either in a sitting or lying position
- The foot and toe measurements are taken without stretching the measuring tape and the toe measurements in particular must be taken with millimeter accuracy
- Toe No 5 is only seldom affected and is therefore only stockinged when an edema is present there
- To be able to measure the foot's length so that the
- measurements are compatible with the circumference measurements, you should mark individual measuring points with a suitable pen, e.g. BELSANA kohl pencil
- Use for measuring toe circumferences the toe measuring tape, which can be simply placed around the single toes

Please see measuring points stipulated in the table.



Measuring point	Circumference (= c)	Length (=I)
(A)	cA ₁ base of instep	IA-A1 inner and outer side of the foot from the metatarsophalangeal joint of the toes to the base of the instep
(A)	circumference from the metatarsophalangeal joint of the CA small toe across the metatarsophalangeal joint of the big toe	
\bigcirc	cX base of the toe at the metatarsophalangeal joint of the toes	IX–Z the length on the inner side of the toe to the centre of the nail. Closed toes are measured to the tip of the toe.
(Z)	cZ centre of the nail	


Taking measurements for flatknitted compression sleeves

Correct and exact measurement is the crucial requirement

• It should be measured in the morning or immediately if a compression arm sleeve is to have a good fit and provide the right medical pressure flow.

- Arms must be free of clothing and free of edema
- Measurements of the arm should be taken on the inside of the arm hanging downwards (the elbow automatically being included in the process)
- upon the removal of a compression bandage
- The circumference should be measured first (mark measuring points with a suitable pen, e.g. BELSANA kohl pencil)

Please see measuring points stipulated in the table. (all measurements in cm)



Measuring point	Circumference (= c)	Length (= I)
$\overline{\mathbb{H}}$	cH from middle of collarbone of one arm to below armpit of the other arm	IG–H from circumference cG to middle of collarbone
G	cG arm below armpit	IC-G from wrist to circumference cG
F	cF widest part of upper arm	IC-F from wrist to circumference cF
E	cE elbow	IC–E from wrist to elbow cE
D	cD widest part of lower arm	IC-D from wrist to circumference cD
(C ₁)	cC ₁ middle of lower arm	IC-C ₁ from wrist to circumference cC1
(C)	cC wrist	IA–C from circumference cA to wrist cC
B	cB hand at height of base of the thumb (measure straight)	IA-B from circumference cA to circumference cB
A	cA hand at height of basal joint of the little finger (measure straight)	


Taking measurements for flatknitted compression gloves

Correct and exact measurement is the crucial requirement if a compression glove is to have a good fit and provide the right medical pressure flow.

- Hand and fingers must be free of edema
- Measurements of the hand should be taken on the inside of the palm
- It should be measured in the morning or immediately upon the removal of a compression bandage

• The circumference should be measured first (mark measuring points with a suitable pen, e.g. BELSANA kohl pencil)

Please see measuring points stipulated in the table.



Measuring point	Circumference (= c)	Length (= I)
W	cW finger at beginning of nail bed	
Y	cY finger at middle joint	IY–W from circumference cY to circumference cW
X	cX finger at basal joint	IX-Y from circumference cX to circumference cY
A	cA palm at height of base of the little finger (measure straight)	
B	cB palm at height of base of the thumb (measure straight)	IA-B from circumference cA to circumference cB
(C)	cC wrist	IA-C from circumference cA to circumference cC
(1)	cC ₁ lower arm at required end point of the glove (depends on individual customer request)	$IA-C_1$ from circumference cA to circumference cC ₁

Range of products (round-knitted)

			SANA osoft		SANA vere		SANA ifortis				SANA v. cottor
		standard I–V	made-to- measure	standard I–V	made-to- measure	standard I–VI	made-to- measure	standard I–VI	made-to- measure	standard I–VI	made-to
	1 18–21 mmHg	•	•			•	•	•	•	•	
Ccl.	2 23–32 mmHg	•	•	•	•	•	•	•	•	•	•
	3 34–46 mmHg									•	•
	A-D	•	•	•	•	•	•	•	•	•	•
	A–F		•		•		•		•		•
	A-G	•	•	•	•	•	•	•	•	•	•
	A-G/T					•	•	• *	• *	•	•
Design	A–T	•	•	•	•	•	•	•	•	•	•
3	A–T/S	•	•	•	•						
	A-T/U	*	•	•	•	•	•		• *		•
	A-T/H					• *	• *		• *		•
	A-T/E						• *		• *		•
		_	_	_	_	_		_		_	
	Silicone	•	•	•	•	•	•	•	•	•	•
	Silicone "fleur"	•	•	•	•	•	•	•	•		•
Grip top	Lace	•	•	•	•			•	•		
	Rubber		•		•		•		•		•
	Effect					•	•				
	long leg (A-D, A-G, A-T)	•		•		•		•		•	
Form	short leg (A–D, A–G, A–T)	•		•		•		•		•	
roiiii	wide fitting (A–D)							•			
	wide fitting (A-G)	•		•		•		•			
	compressive panty p.		•		•		• *		• *		•
	reinforced crotch		•		•		•		•		•
Panty part	comfort gusset		•		•				•	•	•
	vertical slit for men						• *		•		•
	open panty part						•				
	closed toe (one length)	•	•	•	•		•	•	•	•	•
	open toe (one length)	•	•	•	•		•	•	•	•	•
	short and long					•					
Foot	foot closed short and long					•					
	foot open asymm. closed toe		•		•		•		•		
	variable soft toe		-				•		-		

Colours of stockings on the following product pages.

BELSANA grip tops

Lace grip top with silicone (6 cm) – modern refinement

- availbale for BELSANA microsoft and BELSANA vivere
- comfortable grip
- attractive lace appearance with trim for the fashionconscious woman
- slight sheen effect
- breathable due to silicone rhomb surface
- recommended for firm connective tissue

Lace grip top with silicone (6 cm) – feminine elegance

- available for BELSANA classic
- comfortable grip
- no doubling over on account of reinforced borders
- breathable due to striped surface
- recommended for firm connective tissue

Silicone grip top (3 or 5 cm) – classic and discreet

- available for all qualities
- secure grip and shape stability
- discreet, opaque design
- high breathability due to air permeable nub structure
- developed for sensitive skin
- recommended for soft connective tissue

Silicone grip top with motif "fleur" (5 cm) – classic elegance

- available for all qualities
- combines all the advantages of the silicone grip top with a modern look
- opaque design with floral motif

Rubber grip top (3 or 5 cm) – without silicone

- available for all qualities in made-to-measure
- extra secure grip due to anti-slip thread
- no doubling over due to reinforced borders
- discreet striped design
- very high breathability
- ideal in cases of silicone intolerance
- recommended for very soft connective tissue

Effect grip top with silicone (6 cm) – modern striped design

- available for BELSANA comfortis
- comfortable grip
- no doubling over on account of wide striped surface
- breathable due to air channels
- attractive, semi-transparent striped design
- recommended for soft connective tissue









only Ccl. 2



Medical compression stockings:

The rubber stocking? A relict of the past!

perience. Why? Because today's medical compression stockings in no way resemble the thick rubber stockings Wearing medical stockings has never been so pleasurable. of the past. We use highly modern materials and pro-

Compression therapy can be a surprisingly pleasant ex- cesses and offer a wide selection of colours and other fashion elements, such as the elegant lace grip top.

BELSANA medical compression stocking qualities

BELSANA medical compression stockings are available in round-knitted and flat-knitted qualities, each available in various compression classes. A general rule for the choice of a medical compression stocking quality is: the stronger the indication, the more a non-elastic quality is needed.

Roundknitted stockings



BELSANA microsoft supple, fine elegance

Quality

with microfibre Ccl. 1 and 2 Standard sizes I-V and made-to-measure

Especially for

- ideal breathability • fine and transparent quality
- elastic and soft texture modern colours and extras

• permanent skin care with aloe

vera and vitamin E

appearance

comfort

• fine and transparent visual

- Women
- Patients with mild phlebological indications
- Patients with strong connective tissue • Pregnant women
- Fashion-conscious patients



BELSANA vivere with biological aloe vera and vitamin E Ccl. 2

Standard sizes I-V and made-to-measure

BELSANA comfortis soft touch, highly effective Ccl. 1 and 2 Standard sizes I-VI

and made-to-measure

BELSANA classic

 comfort zones ultra-soft material

 maximum wearing comfort • durable yarns

• tested intensively by patients

proven quality for her and for him Ccl. 1 and 2 Standard sizes I-VI soft and fine

opaque appearance

• robust and sturdy quality • large selection of designs

- Patients with dry and sensitive skin • Patients with mild phlebological indications
- Patients with strong connective tissue
- Pregnant women
- soft texture and highest wearing Patients curious about innovative products
 - Patients who are looking for modern and individual solutions
 - Patients with more severe phlebological
 - Patients with slightly soft connective tissue
 - Patients who sit most of the time
 - Patients with more severe phlebological
 - Patients with soft connective tissue
 - Patients who value traditional products



BELSANA classic with cotton natural on the skin

and made-to-measure

Ccl. 1, 2 and 3 Standard sizes I-VI and made-to-measure

- strong and robust quality
- skin-friendly with high cotton content directly on the skin
- great medical efficacy through perfect fit
- Patients who need a high level of compression
- Patients with more severe phlebological indications
- Patients with mild lymphological indications • Patients with soft connective tissue
- Allergy sufferers and lovers of natural fibre

Flatknitted stockings



Quality **BELSANA** impuls

climate regulating pleasant to wear with Merino Ccl. 1, 2, 3 and 4 made-to-measure

Features

- natural high performance fibres from Merino sheep
 - cooling in summer, warming in winter
 - skin-friendly and does not
 - · excellent dimensional stability
- - highly effective working pressure
- **Especially for**
- Patients with mild to severe lipedema and lymphedema • Patients who suffer from burns and scars
- Patients who need individual solutions



BELSANA classic extra

with cotton Ccl. 1, 2 and 3

- verv robust and extra sturdy
- skin-friendly with high cotton content directly on the skin
- perfect fit and high level of pressure
- Patients with mild to severe lipedema and
- Patients who suffer from burns and scars
- · Patients who need individual solutions

made-to-measure



product features

- breathable due to microfibre
- finest fabric envelops the legs
- covers the skin with a feeling of softness
- delicate transparent look
- attractive sheen
- modern colours

Especially for

- Women
- Patients with mild phlebological indications
- Patients with strong connective tissue
- Pregnant women
- Fashion-conscious patients

Material

42 % Polyamide

42 % Elastane

16 % Polyamide (microfibre)

Assortment

Design




stockings



stockings













Tights with patterned panty part

Tights with extra wide panty part

Ccl. 1 and 2

Sizing size I–V

short and long leg length open toe: one foot length closed toe: short and long foot wide fitting (AG with grip top) → see sizing chart p.19

Made-to-measure available

Toe closed, open, asymmetrical closed (made-to-measure)

Grip tops without (elastic cuff), lace (6 cm), silicone (3 cm) (made-to-measure)

silicone (5 cm), silicone with motif "fleur" (5 cm), rubber (3+5 cm) (made-to-measure)

→ see grip top details p.15

normal (A-T), patterned (A-T/S), Panty part

extra wide (A-T/U), compressive (A-T/K) (made-to-measure)

Gusset

Colours

18



additionally: four seasonal trend colours**

* 14 days delivery time

** 14 days delivery time, cannot be returned or replaced

Sizes

75–100	80-110	85-120	95–130	100-135	c ①/U
70-90	75–105	80-110	90-120	95–125	c T
90-115	95-125	100-140	105-145	115–155	c H)/U
90-105	95-115	100-130	105-135	115–145	C (H) A-T Crotch
46-67	49-71	52-74	56-79	58-84	c G for A-T short 69–78 cm long 78,1–88 cm
V	vide fitting only for	thigh length stockin	gs (A-G) with grip to	р	c G for A-G IG
normal wide 46–56 56,1–67	normal wide 49–60 60,1–71	normal wide 52–63 63,1–74	normal wide 56–67 67,1–79	normal wide 58–72 72,1–84	C ((a) for A-G (c) short 65–74 cm long 74,1–84 cm
39-53	41-56	45-60	46-62	47-64	c (F)
31-40	35-43	37-47	39-49	41-52	C E A-D Shaft
28-35	30-38	33-41	35-44	37-47	c D
29–36	32-40	34-43	37-45	39-48	c (C)
23–29	25-32	27-34	30-37	33-40	c Bl
19–20	20,1–22	22,1–24,5	24,6–27	27,1–30	c B
I	II	III	IV	V	w at — with arin ton

w. gt. = with grip top

closed toe open toe

foot length: 15-17,5 cm short foot: 22,5-25,9 cm / shoe size 36-40

long foot: 26-29,5 cm / shoe size 41-45

Product range

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Tights	Tights with patterned panty part	Tights with extra wide panty part
Ccl. 1	Standard Item-No.	218101	-	218104	218141	218162	-
CCI. I	Made-to-measure Item-No.	228101	228103	228104	228141	228162	228143
Cal 3	Standard Item-No.	218201	-	218204	218241	218262	218243
Ccl. 2	Made-to-measure Item-No.	228201	228203	228204	228241	228262	228243

Extras and special designs

	Design	Sizing	Description
grip top A-D, A-F, A-G			
	A-D	mtm	silicone (3+5 cm), rubber (3+5 cm)
grip top	A-F, A-G	standard	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm)
	A-r, A-G	mtm	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm), rubber (5 cm)
panty part			
compressive	A-T, A-T/U	mtm	constant knitted compression (with comfort gusset cH > 138cm)
gusset			
comfort gusset	A-T, A-T/U	mtm	twin seam at the back
reinforced crotch	A-T, A-T/U	mtm	reinforced material inside thigh
toe			
open toe	all	standard/mtm	-
closed soft toe	all	standard/mtm	for standard: short and long foot selectable (see sizing chart)
asymmetrical closed toe	all	mtm	recommended for patients with a long big toe

mtm = made to measure



Product features

- yarn is microencapsulated with 100 % pure aloe vera and vitamin E
- moisturizes the legs consistently even after 100 washing cycles
- leaves a feeling of freshness
- regulates skin temperature effectively
- dry and itchy skin is noticeably improved

Especially for


- Patients with dry and sensitive skin
- Patients with mild phlebological indications
- Patients with strong connective tissue
- Pregnant women
- Patients curious about innovative products

Material

- 42 % Elastane
- 41 % Polyamide (with aloe vera & vitamin E)
- 17 % Polyamide (microfibre)

Assortment

Design




stockings



Half thigh length stockings



Thigh length stockings



Tights with patterned





panty part

Tights with extra wide panty part

Ccl.

Sizing size I–V

short and long leg length open toe: one foot length closed toe: short and long foot wide fitting (AG with grip top) → see sizing chart p.21

Made-to-measure available

Toe closed, open, asymmetrical closed (made-to-measure)

Grip tops without (elastic cuff), lace (6 cm), silicone (3 cm) (made-to-measure)

silicone (5 cm), silicone with motif "fleur" (5 cm), rubber (3+5 cm) (made-to-measure)

→ see grip top details p.15

normal (A-T), patterned (A-T/S), Panty part

extra wide (A-T/U), compressive (A-T/K) (made-to-measure)

Gusset comfort gusset (twin seam at the back) (made-to-measure), reinforced crotch (made-to-measure)

Colours

20

















Sizes

75–100	80-110	85-120	95–130	100-135	c ①/U
70-90	75–105	80-110	90-120	95–125	c T
90-115	95-125	100-140	105-145	115–155	c H/U
90-105	95-115	100-130	105-135	115–145	C H A-T Crotch
46-67	49-71	52-74	56-79	58-84	c for A-T short 69–78 cm long 78,1–88 cm
	wide fitting only for	thigh length stocking	gs (A-G) with grip to	<u>p</u>	c G for A-G IG
normal wide 46–56 56,1–67	normal wide 49–60 60,1–71	normal wide 52–63 63,1–74	normal wide 56–67 67,1–79	normal wide 58–72 72,1–84	c G for A-G is short 65–74 cm long 74,1–84 cm
39–53	41–56	45-60	46-62	47-64	c (F)/
31-40	35-43	37-47	39-49	41-52	C E ID A-D Shaft
28-35	30-38	33-41	35-44	37-47	c D short 34–38 cm long 38,1–43 cm
29–36	32-40	34-43	37-45	39-48	c (C)
23-29	25-32	27-34	30-37	33-40	c (B1)
19–20	20,1–22	22,1–24,5	24,6–27	27,1–30	c B
I	II	III	IV	V	w at = with arin ton

w. gt. = with grip top

open toe closed toe

foot length: 15-17,5 cm short foot: 22,5-25,9 cm / shoe size 36-40

long foot: 26-29,5 cm / shoe size 41-45

Product range

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Tights	Tights with patterned panty part	Tights with extra wide panty part
Cal 2	Standard Item-No.	218001	-	218004	218041	218062	218043
Ccl. 2	Made-to-measure Item-No.	228001	228003	228004	228041	228062	228043

Extras and special designs

mtm = made to measure

	Design	Sizing	Description
grip top A-D, A-F, A-G			
	A-D	mtm	silicone (3+5 cm), rubber (3+5 cm)
grip top	A-F, A-G	standard	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm)
	A-r, A-G	mtm	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm), rubber (5 cm)
panty part			
compressive	A-T, A-T/U	mtm	constant knitted compression (with comfort gusset cH > 138cm)
gusset			
comfort gusset	A-T, A-T/U	mtm	twin seam at the back
reinforced crotch	A-T, A-T/U	mtm	reinforced material inside thigh
toe			
open toe	all	standard/mtm	-
closed soft toe	all	standard/mtm	for standard: short and long foot selectable (see sizing chart)
asymmetrical closed toe	all	mtm	recommended for patients with a long big toe

21

\ /



Product features

- comfort zones at the knee and foot
- ultra-soft material
- maximum wearing comfort and high efficacy
- extremely durable yarns
- opaque appearance
- wide variety of colours
- tested intensively by patients

Especially for


- Patients who are looking for modern and individual solutions
- Patients with more severe phlebological indications
- Patients with slightly soft connective tissue
- Patients who sit most of the time

Material

67 % Polyamide 33 % Elastane

Assortment

Design



Ccl. 1 and 2

Sizing size I–VI

short and long leg length open toe: short and long foot closed toe: short and long foot wide fitting (AG with grip top)

→ see sizing chart p.23

Made-to-measure available

Toe closed, open, asymmetrical closed (made-to-measure), variable (made-to-measure)

Grip tops without (elastic cuff), effect (6 cm), silicone (3+5 cm), silicone with motif "fleur" (5 cm),

rubber (3+5 cm) (made-to-measure), waist belt strap

→ see grip top details p.15

normal (A–T), extra wide (A-T/U), with horizontal or vertical slit (A-T/H) (made-to-mea-Panty part

sure), compressive (A-T/K) (made-to-measure)

reinforced crotch (made-to-measure) Gusset

Colours

22



Sizes

75–100	80-110	85-120	95–130	100-135	105-140	c (T)/U
70–90	75–105	80-110	90-120	95–125	100-130	
90-115	95–125	100-140	105–145	115–155	120–165	c H/U
90-105	95–115	100-130	105–135	115–145	120–155	C (H) A-T Crotch
46-67	49-71	52-74	56-79	58-84	62-89	c G for A-T short 69–78 cm long 78,1–88 cm
	wide fittin	g for thigh length	stockings (A-G) wi	th grip top		c G forA-G
normal wide 46–56 56,1–67	normal wide 49–60 60,1–71	normal wide 52–63 63,1–74	normal wide 56–67 67,1–79	normal wide 58–72 72,1–84	normal wide 62–76 76,1–89	A-G Leg
39–53	41-56	45-60	46-62	47-64	49-68	C F
31-40	35-43	37-47	39-49	41-52	43-56	c E
28-35	30-38	33-41	35-44	37-47	39–51	C D A-D Shaft short 34–38 cm
29-36	32-40	34-43	37-45	39-48	41-52	c C long 38,1–43 cm
23-29	25-32	27-34	30-37	33-40	37-44	c B1
19–20	20,1–22	22,1–24,5	24,6–27	27,1–30	30,1–33	c B
28-32	29-34	30-35	32–36	33-38	34-39	c (Y)
19-23	20-24	21-25	22–26	23-27	24-28	c (A)
1	II	Ш	IV	V	VI	w. gt. = with grip top

open toe closed toe

short foot: 15-18,4 cm short foot: 22,5-25,9 cm / shoe size 36-40 long foot: 18,5-22 cm long foot: 26-29,5 cm / shoe size 41-45

Product range

mtm = made to measure

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Thigh length stockings with belt strap	Tights	Tights with extra wide panty part	Tights for men with horizontal slit	Tights with one leg
Cal 1	Standard Item-No.	215101	-	215104	-	215141	215143	-	-
Ccl. 1	Made-to-measure Item-No.	225101	225103	225104	225147	225141	225143	-	-
ددا ۲	Standard Item-No.	215301	-	215304	215347	215341	215343	215346	-
Ccl. 2	Made-to-measure Item-No.	225301	225303	225304	225347	225341	225343	225346	225345

	Design	Sizing	Description
grip top A-D, A-F, A-G			
	A-D	standard	effect (6 cm), silicone (3+5 cm)
avia tan	A-D	mtm	effect (6 cm), silicone (3+5 cm), rubber (3+5 cm)
grip top	A E A G	standard	effect (6 cm), silicone (5 cm), silicone "fleur" (5 cm)
	A–F, A–G	mtm	effect (6 cm), silicone (5 cm), silicone "fleur" (5 cm), rubber (5 cm)
panty part			
compressive	A-T, A-T/U	mtm	constant knitted compression. Available in Ccl. 2
with vertical slit	A-T/H	mtm	vertical slit
with stoma opening	A-T	mtm	location neccessary
open panty part	A-T/U	mtm	open panty part front
gusset			
open crotch	A-T	mtm	open crotch area
reinforced crotch	all panty parts	mtm	reinforced material inside thigh
toe			
open toe	all	standard/mtm	for standard: short and long foot selectable (see sizing chart)
closed soft toe	all	standard/mtm	for standard: short and long foot selectable (see sizing chart)
asymmetrical soft toe	all	mtm	recommended for patients with a long big toe
variable soft toe	all	mtm	length individual selectable (5–10 cm measuring details see p. 10). Recommended for Hallux Valgus



Product features

- extra flexible material
- multifilament fibres make the stocking robust and breathable
- large selection of stocking types to choose from
- semi-opaque appearance

Especially for


- Patients with more severe phlebological indications
- Patients with soft connective tissue
- Patients who value traditional products

Material

74 % Polyamide 26 % Elastane

Assortment

Design



Ccl. 1 and 2

Sizing size I–VI

short and long leg length open toe: one foot length closed toe: short and long foot

wide fitting (AD without grip top, AG with grip top)

→ see sizing chart p.25

Made-to-measure available

Toe closed, open, asymmetrical closed (made-to-measure)

Grip tops without (elastic cuff), lace (6 cm), silicone (3+5 cm), silicone with motif "fleur" (5 cm),

rubber (3+5 cm) (made-to-measure), waist belt strap

→ see grip top details p.15

Panty part normal (A–T), extra wide (A–T/U) (made-to-measure), compressive (A–T/K) (made-to-measure),

with vertical or horizontal slit (A–T/H) (made-to-measure)

Gusset comfort gusset (twin seam at the back) (made-to-measure), reinforced crotch (made-to-measure)

Colours

24



Sizes

62-84	67–89	74–98	80-102	85–108	90-110	c (T)	
70–95	77–100	84–116	90–127	98–137	105–148	c H	IT
42-55	44–58	48–62	50–67	53–73	56–80	C G for A–T and A–G w/o gt	IK
wide fitting fo	r thigh length sto	ckings (A–G), com	nbined with norm	al circumferences	for cB1, cC, cD	c G for A-G wide	A–T Crotch short 65–74 cm
normal wide 44–53 51–60	normal wide 46–56 54–63	normal wide 50–60 58–68	normal wide 53–64 63–74	normal wide 57–70 67–79	normal wide 61–78 71–86	C G for A–G normal w.gt.	long 74,1–84 cm
normal wide 38–48 47–55	normal wide 41–51 49–58	normal wide 45–54 53–62	normal wide 46–57 56–66	normal wide 47–60 58–69	normal wide 49–62 59–71	C F for A-G wide C F IF	A-G Leg short 65–74 cm long 74,1–84 cm
29–39	31–41	33-43	36–46	38–49	40-52	C (E)	A-G/T Leg long 74,1–84 cm
	wide fitting or	nly for knee high s	tockings (A–D) w	ithout grip top		c D for A–D wide	ID
normal wide 26–33 29–37	normal wide 28–36 32–40	normal wide 30–39 34–44	normal wide 33–42 39–48	normal wide 35–45 41–51	normal wide 37–48 43–54	c D for A-D IC	A–D Shaft short 34–38 cm long 38,1–43 cm
normal wide 28–34,5 32–37	normal wide 30–37 35–41	normal wide 32–40 38–45	normal wide 35–43 40–49	normal wide 37–46 42–52	normal wide 39–49 44–55	c for A-D IB1	LA
normal wide 23–27 25–29	normal wide 25–29 27–31	normal wide 27–31 29–33	normal wide 29–33 32–36	normal wide 31–36 34–39	normal wide 35–40 38–43	c B for A-D romal	IA
19–20,5	20,6-22,5	22,6-24,5	24,6-27	27,1–30	30,1–33	C Holling , ,	= with grip top
26-30	27-32	29-34	30-37	32-39	33-41	c (Y)//,// w/o gr	t = without grip top
19–20,5	20,5-22,5	22,5-24,5	24,5–27	27-30	30-33	c (A)/	
I	II	III	IV	V	VI		

open toe

closed toe

Important information as to the wide fitting option

foot length: 15–17,5 cm short foot: 22,5–25,9 cm / shoe size 36–40 long foot: 26–29,5 cm / shoe size 41–45

Wide circumference A-G cannot be combined with wide circumference A-D, for wide circumference A-G normal measurements apply in the case of cB1, cC and cD!

Product range

mtm = made to measure

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Thigh length stockings with belt strap	Tights	Tights with extra wide panty part	Tights for men with horizontal slit	Tights with one leg
Col 1	Standard Item-No.	215201	-	215204	-	215241	-	-	-
Ccl. 1	Made-to-measure Item-No.	225201	225203	225204	-	225241	-	-	-
Calla	Standard Item-No.	215401	-	215404	215447	215441	-	-	-
Ccl. 2	Made-to-measure Item-No.	225401	225403	225404	225447	225441	225443	225446	225445

	Design	Sizing	Description
grip top A-D, A-F, A-G			
	A-D	standard	silicone (3+5 cm)
arin ton	N-D	mtm	silicone (3+5 cm), rubber (3+5 cm)
grip top	A F A C	standard	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm)
	A–F, A–G	mtm	lace (6 cm), silicone (5 cm), silicone "fleur" (5 cm), rubber(5 cm)
panty part			
compressive	A–T, A–T/U, A–T/H	mtm	constant knitted compression. Available in Ccl. 2
with vertical slit	A-T/H	mtm	vertical slit
with stoma opening	A-T	mtm	location neccessary
gusset			
comfort gusset	A-T, A-T/U	mtm	twin seam at the back
open crotch	A-T	mtm	open crotch area
reinforced crotch	all panty parts	mtm	reinforced material inside thigh
toe			
open toe	all	standard/mtm	-
closed soft toe	all	standard/mtm	for standard: short and long foot selectable (see sizing chart)
asymmetrical closed toe	all	mtm	knit to form. Recommended for long bunion

27



Product features

- high cotton content directly on the skin
- extremely robust
- shape-retaining material exerts deep pressure on the legs
- natural cooling effect due to cotton
- perfect for allergy sufferers

Especially for

- Patients who need a high level of compression
- Patients with more severe phlebological
- Patients with mild lymphological indications
- Patients with soft connective tissue
- Allergy sufferers and lovers of natural fibre

Material

46 % Polyamide

31 % Elastane

23 % Cotton

Assortment

Design



Ccl. 1, 2 and 3

Sizing size I–VI

short and long leg length open toe: one foot length closed toe: short and long foot

→ see sizing chart p.27

Made-to-measure available

Grip tops without (elastic cuff), silicone (3+5 cm), silicone with motif "fleur" (5 cm) (made-to-measure),

rubber (3+5 cm) (made-to-measure), waist belt strap

→ see grip top details p.15

normal (A-T), extra wide (A-T/U) (made-to-measure), compressive (A-T/K) Panty part

(made-to-measure), with vertical or horizontal slit (A–T/H) (made-to-measure)

comfort gusset (twin seam at the back) (Standard), reinforced crotch (made-to-measure) Gusset

Colours

26

Toe



mode bright mode









Sizes

26–30 19–20,5	27–32 20,5–22,5	29–34 22,5–24,5	30–37 24,5–27	32–39 27–30	33–41 30–33	c (Y) IB IA
19–20,5	20,6–22,5	22,6–24,5	24,6–27	27,1–30	30,1–33	A-D Shaft short 34–38 cm long 38,1–43 c
23-27	25-29	27-31	29-33	31–36	35-40	C (B1)
28-34,5	30-37	32-40	35-43	37-46	39-49	C (C) IE long 74,1–84 c
26-33	28-36	30-39	33-42	35-45	37-48	c D long 74,1-84 c
29-39	31–41	33-43	36–46	38-49	40-52	C (E) A-G Leg short 65-74 cn
38-48	41-51	45-54	46-57	47-60	49-62	_ c (F)\(\)
44-53	46-56	50-60	53-64	57-70	61-78	c G w. gt short 65-74 cn long 74,1-84 cr
42-55	44-58	48-62	50-67	53-73	56-80	C G A-G w/o. gt. K
70-95	77-100	84-116	90-127	98-137	105-148	c (H)
62-84	67-89	74-98	80-102	85-108	90-110	_ c (T)

closed toe open toe

foot length: 15-17,5 cm short foot: 22,5-25,9 cm / shoe size 36-40

long foot: 26-29,5 cm / shoe size 41-45

w.gt. = with grip top w/o gt = without grip top

Product range

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Thigh length stockings with belt strap	Tights	Tights with extra wide panty part	Tights for men with horizontal slit	Tights with one leg right / left
Ccl. 1*	Standard Item-No.	213301	-	213304	213347	213341	-	-	-
CCI. I	Made-to-measure Item-No.	-	-	-	-	-	-	-	-
Ccl. 2	Standard Item-No.	213401	-	213404	213447	213441	-	-	-
CCI. Z	Made-to-measure Item-No.	223401	223403	223404	223447	223441	223443	223446	223445
C-I 2	Standard Item-No.	213501	-	213504	213547	213541	-	-	-
Ccl. 3	Made-to-measure Item-No.	223501	223503	223504	223547	223541	223543	223546	223545
							* Ccl. 1 is avail	able with open too	e in colour mode

	Design	Sizing	Description
grip top A-D, A-F, A-G			
	A-D	mtm	silicone (3+5 cm), rubber (3+5 cm)
grip top	A–F, A–G	standard	silicone (5 cm)
	Λ-I, Λ-G	mtm	silicone (5 cm), silicone "fleur" (5 cm), rubber (5 cm)
panty part			
compressive	A–T, A–T/U, A–T/H	mtm	constant knitted compression. Available in Ccl. 2 and 3
with vertical slit	A-T/H	mtm	vertical slit
with stoma opening	A-T	mtm	location neccessary
gusset			
open crotch	A-T	mtm	open crotch area
reinforced crotch	all panty parts	mtm	reinforced material inside thigh
toe			
open toe	all	standard/mtm	-
closed toe	all	standard/mtm	standard: short and long foot selectable (see sizing chart)

nach Maß



climate regulating, pleasant to wear

Flatknitted compression stockings with the natural properties of Merino wool

Product features

- natural high performance fibres from Merino sheep
- excellent dimensional stability
- breathable and moisture regulating
- cooling in summer, warming in winter
- robust, durable and odour neutral
- skin-friendly and does not scratch
- pleasantly soft due to very fine
- highly effective working pressure

Especially for


- Patients with lipedema and lymphedema
- Patients who suffer from burns and scars

Material

64 % Polyamide (Microfibre), 23 % Elastane, 13 % Merino wool

Assortment

Design



Ccl. 1, 2, 3 and 4

Made-to-measure only made-to-measure

Toe closed, open, asymmetrical foot closed, asymmetrical foot open

Heel

without (elastic cuff), silicone (3 + 5 cm), silicone with motif "fleur" (5 cm), **Grip tops** rubber (3+5 cm), grip top pieces (3 + 5 cm), Flex-Band (5 cm) (A-D only)

lateral elevation, anatomical knee, stump, upper leg elevation, with stirrup Leg

normal (A-T), extra wide (A-T/U), vertical or horizontal slit (A-T/H), lateral elevation (slipform), Panty part

hip band, elastic band (3 cm)

comfort gusset (twin seam at the back), breathable gusset, open crotch, reinforced crotch Gusset

Fastener zip fastener short and long

Pads lymph pads DIN A4, pressure pad, with or without pocket, fabric lining

Colours

















Deco bands

28





Product range

		Knee high stockings	Half thigh length stockings	Thigh length stockings	Thigh length stockings with belt strap	Tights	Tights with extra wide panty part	Tights for men with horizontal slit
Ccl. 1	Made-to-measure Item-No.	222501	222503	222504	222547	222541	222543	222546
Ccl. 2	Made-to-measure Item-No.	224501	224503	224504	224547	224541	224543	224546
Ccl. 3	Made-to-measure Item-No.	226501	226503	226504	226547	226541	226543	226546
Ccl. 4	Made-to-measure Item-No.	227501	227503	227504	227547	227541	227543	227546
		Tights with	Leg section	Leggings	Capri pants	Bermuda	Cyclist	
		one leg	Leg section (C–G)	33 3		shorts	trousers	
Ccl. 1	Made-to-measure	one leg	•	222550	222552	shorts 222554	•	
Ccl. 1		one leg	Leg section (C-G)				trousers	
	Item-No. Made-to-measure	one leg right / left 222545	Leg section (C-G) 222518	222550	222552	222554	trousers 222555	

	Design	Description
grip top A-D, A-F, A-G, C-C		
arin tons	A-D	silicone (3+5 cm), silicone "fleur" (5 cm), rubber (3+5 cm), Flex-Band (5 cm)
grip tops	A–F, A–G and C–G	silicone (5 cm), silicone "fleur" (5 cm), rubber (5 cm)
grip top pieces	all	silicone (3 + 5 cm). Variable positioning or ¾ cG circumference
toe		
open	all	
closed	all	
asymmetrical foot closed	all	anatomical foot-shaped design structured for better lymphatic flow
asymmetrical foot open	all	anatomical foot-shaped design structured for better lymphatic flow
heel		
90°-heel	all	shape of foot at 90° angle
leg		
lateral elevation	A-D, A-F, A-G and C-G	elevation on the side of the leg
upper leg elevation	A-G	elevation on the front of the leg
anatomical knee	A–F and higher	anatomical angularity of the knee area
stump	all	open or closed. For patients with leg amputation
with stirrup	B-T	rubber band to fix the leggings
panty part		
slip design	all panty parts	elevation on the side of the leg
vertical slit	A-T/H	vertical slit for men
hip band	all panty parts	knitted on waist trim
elastic band (3 cm)	A–T	soft waist band with snug fit. Prevents unpleasant folding down of top
gusset		
comfort gusset	A-T, A-T/U	twin seam at the back
breathable gusset	all panty parts	gusset with permeable perforation
open crotch	A–T	gusset with recess
reinforced crotch	all panty parts	reinforced material inside thigh
zip-fastener (incl. hooks)		
zip fastener short	A-D	indicate length and positioning. Zip on inside, outside and in the seam possible
zip fastener long	A-F, A-G	indicate length and positioning. Zip on inside, outside and in the seam possible
pads		
lymph pad	all	pad with profile made of foam material for improved lymphatic flow and selective increase in pressure. Available in DIN A4 for trimming to size or pre-trimmed to size
pressure pad	all	Anatomically moulded silicone pad for optimum pressure distribution on the inside and outside ankle
lymph pad or pressure pad sewn in	all	lymph pad or pressure pad are securely integrated in the stocking
lymph pad or pressure pad pocket sewn in	all	lymph pad or pressure pad can be, if required, removed from pocket
lining fabric sewn in	A–F and higher	material sewn in on inside of stocking for enhanced comfort
asseccories		
deco bands	A-D, A-G, A-G/T, A-T, A-T/U	coloured seams in croco, silver and bronze

classic extra

individual care

Flatknitted compression stockings with comfortable cotton for arm, hand and toe

Product features

- skinfriendly and robust thanks to high cotton content
- strong and shape-retaining material exerts effective deep
- individual treatment
- extra flat seams

Especially for

- Patients with mild to severe lipedema and lymphedema
- Patients who suffer from burns

Material

40 % Cotton

28 % Polyester

23 % Elastane

9 % Polyamide

Assortment

Design





Forearm-hand Upper arm-hand Shoulder-hand without thumb



(seam inside) (seam inside)



without thumb without thumb (seam inside)



Upper arm forearm



forearm

thumb

Glove without



Glove with

thumb



Glove with (open

or closed) fingers









Foot and toe part (open or closed toes)

Ccl. 1, 2, 3

Made-to-measure only made-to-measure

Fingers/Toes

Grip tops without (elastic cuff), silicone (3 + 5 cm), rubber (3 + 5 cm)

closed, open, with thumb, without thumb

Holder sleeve holding strap (A-H, C-H), bra holder (A-H, C-H), silicone grip top pieces (3 + 5 cm)

lateral elevation, stump stocking Arm

Pads lymph pads DIN A4, pressure pad sewn in

Colours

Fastener





zip fastener short and long

Product range

		Elbow section (D–F)	Forearm – hand without thumb (A–D) (seam inside)*	Upper arm – hand without thumb (A–G) (seam inside)*	Shoulder – hand without thumb (A–H) (seam inside)*	Upper arm – forearm (C – G)	Shoulder- forearm (C-H)
Ccl. 1	Made-to-measure Item-No.	222471	222473	222475	222474	222477	222476
Ccl. 2	Made-to-measure Item-No.	224471	224473	224475	224474	224477	224476
Ccl. 3	Made-to-measure Item-No.	226471	226473	226475	226474	226477	226476

^{*} two parted treatment recommended, one parted treatment only available in skin

		Glove without thumb	Glove with thumb	Glove with fngers (open or closed)	Finger (open or closed)	Foot and toe part (open or closed)
Ccl. 1	Made-to-measure Item-No.	222470	222480	222481	222483	222494
Ccl. 2	Made-to-measure Item-No.	224470	224480	224481	224483	224494
Ccl. 3	Made-to-measure Item-No.	226470	226480	226481	226483	226494

	Design	Description
grip top A-D, A-G, C-G, D-F		
grip top	A-D, A-G, C-G, D-F	silicone (3+5 cm), rubber (3+5 cm)
grip top pieces	all	silicone (3+5 cm), individual positioning
holder for A-H & C-H		
sleeve holding strap	A-H, C-H	strap is fixed diagonally across the chest
bra holder	A-H, C-H	to be hooked into the bra strap
arm		
lateral elevation	A-D, A-G, C-G	elevation outer arm
stump	all	for patients with arm amputation
zip-fastener		
zip fastener short	all	indicate length and positioning
zip fastener long	all	indicate length and positioning
pads		
lymph pad DIN A4	all	pad with profile made of foam material for improved lymphatic flow and selective increase in pressure
pad sewn in	gloves	silicone pad for gloves

Accessories and Care

BELSANA magnide® – donning aid for medical compression stockings

- for open and closed toe
- two in one: magnide® on/off technology for fast and easy use
- handy, practical and space-saving perfect for on the way
- protects fabric
- reduction of physical strain for both caregiver and stocking wearer
- increase in independence of the stocking wearer
- increase in therapy compliance
- available in 2 sizes

Check the BELSANA youtube channel to view donning videos

BELSANA Angel – donning aid for medical compression stockings

- ideal when personal mobility is restricted
- stable frame with two moveable handles
- perfect for every day use at home
- easy to use, easy to clean
- suitable for all designs with closed and open toe

Check the BELSANA youtube channel to view donning videos

BELSANA grip-Star special gloves with grip

- facilitates the putting-on and taking-off of medical compression stockings
- secure grip thanks to profile on the palm side of the glove
- protects the knitted fabric against accidental
- latex coat without direct skin contact, can therefore be used by allergy sufferers
- Material: 55 % Polyester, 45 % Latex

BELSANA hafti - skin adhesive for medical compression stockings

- can be used as a grip-top alternative
- gentle grip on the skin
- ph-neutral and skinfriendly
- with natural ingredient camomille
- water-soluble

Content	1 pi	ece
Item-No.	29980508	29980512
Shoe size	≤ 40	41-45
Size	medium	large
		I



Please note: BELSANA magnide® contains several magnets. Please keep it at a minimum distance of 30 cm from pacemakers and defibrillators.

	Angel small cC<44cm	Angel big cC ≥ 44 cm
Item-No.	299809	299808
Content	1 piece	1 piece



Content			1 pair		14	
Item-No.			299818		1	-
Size	extra-small	small	medium	large	extra-larg	ge

Please note: the Belsana grip-Star facilitates use and helps handle the BELSANA donning aids better.

Item-No.

Content



BELSANA suspender single / double

- worn across the shoulders or for double
- suitable for compression thigh lengths stockings without grip top
- colour: skin

	single	double
Item-No.	299830	299815
Content	1 piece	





BELSANA panty belt

- with 2 pairs of suspenders on each side
- suitable for compression thigh lengths stockings
- colour: skin

Circumference	80 cm	90 cm	100 cm	110 cm	120 cm	130 cn
Item-No.			299850			
Content			1 piece			



BELSANA sleeve holding strap

• suitable for BELSANA classic extra A–H and C–H sleeve segments

Item-No.	224484
Content	1 piece



BELSANA bra holder

• suitable for BELSANA classic extra A-H and C-H sleeve segment

Item-No.	224485
Content	1 piece



BELSANA blau40 - special detergent for medical compression stockings

- washes gently, maintains fabric and freshens up colours
- skin-friendly and deodorant

Item-No.	2998005
Content	1 bottle à 250 ml







BELSANA washing net

- for gentle washing of medical compression stockings
- size: 300 x 360 mm
- Material: 100 % Polyester

Item-No.	298096
Content	1 piece




32



- protect against venous leg problems
- active support around ankles
- soft stocking top for comfortable hold

Sizes				
	- 1	II	III	I۱
Shoe sizes	36-38	39-41	42-44	45-



BELSANA cotton

with comfortable cotton

50 % Cotton, 40 % Polyamide, 10 % Elastane


BELSANA sommer Edition

53 % Polyamide, 24 % Polyamide (microfibre)

with breathable microfibre

13 % Cotton, 10 % Elastane (Lycra)

Knee high stockings A-D



Knee high stockings A-D

Item-No. 232101







Content

1 pair

Support stockings:

Freedom for your legs!

Swollen feet? Legs heavy as lead? Not only women experiorrect pressure distribution that decreases from the foot ence this problem: men also recognize the symptoms after standing or sitting at work all day. Or after a long journey by plane or by car. Long periods of standing or sitting can trigger blood congestion in the veins, leading to a gradual in a wide variety of colours and designs. enlargement of the vein walls. Support stockings with the

upwards promote the flow of blood and are therefore an ideal preventative measure for both men and women. BELSANA support stockings are fashionable and available

with moisturising aloe vera

42 % Polyamide (with aloe vera), 35 % Polyamide, 15 % Cotton, 8 % Elastane

Knee high stockings A-D











BELSANA vital Edition

Colours









Permanent skin care with 100 % pure, certified organically grown aloe vera (micro-encapsulated). Refreshing effect thanks to high moisture content – soft touch.

traveller po



Compressive knee-high stockings especially for travelling

Product features

- correct, decreasing graduated pressure upwards along the leg
- individual size determined on the basis of two measuring points (see sizing chart)
- discreet ribbed look
- non-constricting comfort cuff
- medically effective
- a preventive measure against travel thrombosis
- offers protection against swollen feet and legs
- recommended by frequent flyers
- unisex for her and for him


Colours



Material 84 % Polyamide,

16 % Elastane

Design



stockings



izes					
	small	medium	large	x-large	
circumference cm (cC)	30–36	34–41	39–47	43–51	©
circumference cm (cB)	19–21	21,1-24	24,1–27	27,1–31	B
	1	II	III	IV	
Shoe sizes	35–38	39–42	43-46	47–50	

Knee high stockings A-D

Content	Item-No.
1 pair	231701

traveller woman



Compressive knee-high and thigh length stockings for travelling especially for women

Product features

Colours

- correct, decreasing graduated pressure upwards along the leg
- light sheen, fine transparency
- sock size based on individual measurements (see sizing chart)
- a preventive measure against travel thrombosis
- offers protection against swollen feet and legs
- recommended by frequent flyers

Material

52 % Polyamide (microfibre), 29 % Elastane, 19 % Polyamide

Design



Knee high stockings



Sizes		
	medium	
circumference cm (cG)	50-63	c (G)
circumference cm (cC)	34–41	c ©
circumference cm (cB)	21,1–25	c B

Length (cm)	standard	long	
IG A–G Leg	66,1–75	75,1–84	
ID A–D Shaft	34,1–38	38,1–43	
Shoe sizes	36–42		

Knee high stockings A-D

Content	Item-No.

Thigh length stockings A-G

1 pair	233804T
Content	Item-No.

70den

Knee high stockings, thigh length stockings and tights with active support

with extra wide

panty part

Product features


- decreasing support strength from the ankle upwards
- perfect fitting
- attractive colours
- very transparent look
- moulded heel
- small flexible soft toe



Material

80 % Polyamide, 20 % Elastane

Design



Knee high stockings A-D

Comfortable soft cuff

1 pair	233101
Content	Item-No.

Sizes

	1	II	III
Shoe sizes	35–37	38–40	41–43

stockings

Thigh length stockings A-G

Attractive lace grip top (7 cm)

1 pair	233104
Content	Item-No.

Sizes

	_	1				-
Dress size	36-38	38-40	40-42	42-44	44-46	48-50

Tights A-T

Modern lace panty part, size VI with comfort gusset

Content	Item-No.
1 piece	233141

Sizes

	0	- 1	II	III	IV	V	VI
Dress size	36–38	38-40	40-42	42-44	44–46	48-50	50-52

Tights with extra wide panty part A-T/U

Comfort panty part, ideal for full-figured persons and for pregnant women, available colours: lava, bronze

Content	Item-No.
1 piece	233143

Sizes

	0	1	II	III	IV	V
Dress size	36–38	38-40	40-42	42-44	44-46	48-50

140den

Knee high stockings, thigh length stockings and tights with high active support

Product features

- decreasing support strength from the ankle upwards
- perfect fitting
- small flexible soft toe
- moulded heel
- attractive colours
- soft feel, fine and slight transparent look


Colours



Material

66 % Polyamide, 34 % Elastane

Design



stockings

stockings

Sizes

Comfortable soft cuff

Knee high stockings A-D

Content	Item-No.
1 pair	233301

	1	II	III
Shoe sizes	35–37	38–40	41–43

Thigh length stockings A-G

Attractive lace grip top (7 cm)

Content	Item-No.
1 pair	233304

Sizes

	0	1	II	III	IV	٧
Dress size	36–38	38–40	40-42	42-44	44–46	48–50

Tights A-T

Modern lace panty part, sizes IV and V with comfort gusset

1 piece	233341
Content	Item-No

	0	1	II	III	IV	٧
Dress size	36-38	38-40	40-42	42-44	44-46	48-50

open toe

Strong support – glossy look!


Product features

- decreasing support strength from the ankle upwards
- support complies with regulations for medical compression stockings
- moulded reciprocated heel and reinforced toes
- size determination according to cB measure

Material

52 % Polyamide (microfibre), 29 % Elastane, 19 % Polyamide

Design





stockings







Tights with extra wide panty part

Colours











stockings



additionally: four seasonal trend colours*

open toe is available in the following colours: siena, pearl, night blue, black (5 days delivery time)

Sizes

Circumference (cm)

	18-	-21	21,1	-25	25,1	-30	C	B	
								_ / / / /	
	29-	-35	34-	-41	37-	-49	С	©	
A-G with grip top	normal 43–52	wide 52–60	normal 50–63	wide 63-74	normal 55–73	wide 73–90	С	(G w.gt. /	D N-D
for tights	43-	-60	50-	-74	55-	-90	С	G for A-T.	∖–G/
extra wide panty part		-120 -100	105- 95-	-130 110	115- 100-	-150 -120	C	H)/U/	G
extra wide panty part		-88 -80	79– 70-		95- 85-		С	①/U	

	Length (cm)	short	standard	long
	IG A-G and A-T leg	58–66	66,1–75	75,1–84
	ID A-D shaft	32–34	34,1–38	38,1–43
ľ				

36-42

w. gt. = with grip top

Knee high stockings A-D

1 pair	233801
Content	Item-No.

Thigh length stockings A-G with lace grip top (7 cm) normal

1 pair	233804
Content	Item-No.

Thigh length stockings A-G with lace grip top (7 cm) wide fitting

Content	Item-No.
1 pair	233804

Tights A-T with modern panty part

Shoe sizes

1 piece	233841
Content	Item-No.

Tights A-T/U

with extra wide panty part

ideal for full-figured persons and pregnant women

Content	Item-No.	
1 piece	233843	

* 14 days delivery time, cannot be returned or replaced



SUPPORT STOCKINGS

PERFORMANCE



PERFORMANCE

True Medical Compression for the beginner and the more advanced sportsperson

compression socks low cut

Product features

- support for ligaments and tendons
- secure hold and stability
- fewer injuries
- functional shock absorption due to high-tech plush sole
- wrinkle-free fit, no friction or blisters
- moisture wicking functional fibres
- silver finish for odour reduction and anti-bacterial

Colours







Socks A-B

Content	Item-No.
1 pair	232705

Material

70 % Polyamide, 20 % Polypropylene, 10 % Elastane

Design





	- 1	II	III	IV
Shoe sizes	36–38	39–41	42-44	45-47

compression socks high

Product features

- with effective compression technology
- anatomically correct compression
- support venous functions and blood circulation
- support for ligaments and tendons
- secure hold and stability
- fewer injuries
- functional shock absorption due to high-tech plush sole
- wrinkle-free fit, no friction or blisters
- moisture wicking functional fibres
- silver finish for odour reduction and anti-bacterial

Colours







Knee high stockings A-D

Content	Item-No.
1 pair	232701

Material

58 % Polyamide, 29 % Polypropylene, 13 % Elastane

Design





Sizes

	1	II	III	IV
Shoe sizes	36–38	39–41	42–44	45–47

PROFESSIONAL



True Medical Compression for the professional and the competitive athlete

compression socks medium cut

Product features

- with professional compression technology
- fewer injuries thanks to stabilizing bandaging function
- ideal for all sports with abrupt movements
- selective padding of inner and outer ball of the foot
- asymmetric anatomically matched sole (r/l)
- less fatigue and faster regeneration
- individual fit: individually measured and fitted (see sizing chart)
- profile retaining and functional shock absorption thanks to high tech yarns in the sole
- specially designed anti-slip honeycomb sole for an optimal hold in the shoes
- toe and heel padding to prevent blisters
- seamless design for less friction
- high breathability through high-tech climate fibres

Colours



Socks A-B.

1 pair	234605
Content	Item-No.

Material

59 % Polyamide, 26 % Polyester, 15 % Elastane

Design







	small	medium	large	x-large	
ircumference cm :B)	19–21	21,1-24	24,1-27	27,1–31	(

	I	Ш	III	IV
Shoe sizes	35–38	39–42	43-46	47–50

PROFESSIONAL



True Medical Compression for the professional and the competitive athlete

compression socks high

Product features

- with professional compression technology
- ideal for all running sports and ball games
- faster warm up and regeneration of the muscles before and after sports
- supports blood circulation and venous function
- individual fit: individually measured and fitted
- size determination according to 2 measuring points
- toe and heel padding
- · forefoot stabilization
- secure hold and stability to reduce risk of injuries
- asymmetric anatomically matched sole (r/l)
- anatomical foot and ankle padding to prevent
- profile retaining and functional shock absorption
- anti-bacterial and odour reducing silver ion yarn
- moisture wicking climate fibre

Knee high stockings A-D

1 pair	234601		
Content	Item-No.		

Material

83 % Polyamide, 17 % Elastane







Design



		medium	•	•	
circumference cm (cC)	30-36	34–41	39–47	43–51	c (
circumference cm (cB)	19–21	21,1–24	24,1–27	27,1–31	C
	- 1	II	III	IV	

compression tube

Product features

- individual graduated compression
- supports blood supply of muscles and tissue
- individual fit: individually measured and fitted (see sizing chart)
- size determination according to 2 measuring
- performance enhancing effect during sports
- less fatigue and faster regeneration
- more efficient muscle capacity thanks to reduced muscle vibration
- ideally worn under training wear, wetsuits or barefoot
- Fine knitted fabric with soft rim


Tube D-B


Content	Item-No.
1 pair	234621

Material

Shoe sizes

62 % Polyamide, 38 % Elastane







COLOURS

39–42 43–46 47–50

Sizes

	small	medium	large	x-large	
circumference cm (cC)	30–36	34–41	39–47	43-51	c (C)
circumference cm (cB)	19–21	21,1–24	24,1–27	27,1–31	CB

more information under www.belsana-sport.com

soft

Special socks for diabetics with antimicrobial silver fibre

Product features

- without support effect
- temperature-regulating
- soft cuffs and external seams, do not exert pressure

Socks A-B.

Content	Item-No.
1 pair	232406

Material

68 % Cotton, 29 % Polyamide, 3 % Silver





iizes				
	1	II	Ш	IV
Shoe sizes	36–38	39–41	42-44	45–47

The silver fibre in the BELSANA diabetic sock has special properties that are ideally suited to help prevent typical foot problems suffered by diabetics.

Hygiene through anti-bacterial effect of the silver

The silver fibre content transforms one's sweat. Bacteria reduction is achieved and, as a result, speed of the reduction process plays a major reproducing and spreading themselves.

Optimum foot climate

Silver is the element with the highest degree of conductivity. This property makes it possible that the foot climate is regulated in a nathere is no unpleasant odour build-up. The utral way. Body heat is stored during the cold directly on the skin and ensure a sustained season and returned to the body. In the case of role, as germs and fungi are prevented from warm outside temperatures body heat is transferred outwards.

Outstanding wearing comfort

points almost totally knitted without seams. The antimicrobial properties of silver develop feeling of freshness. The soft trim of the special diabetic sock also creates a most pleasant feeling of comfort. Constricting rubber threads have been totally dispensed with.

BELSANA soft is in order to prevent pressure



Terms and conditions of sale and delivery

General: All deliveries, services and others shall be performed on the basis of these terms and conditions of trade. Purchaser's terms and conditions that deviate from these, and which have not been explicitly acknowledged in writing by the supplier, shall not be considered binding on the supplier, even if the supplier has not explicitly objected to them. The rights and obligations arising from the business relationships with the supplier may not be passed on to third parties.

Prices

Prices are quoted in Euro. Prices are net prices and quoted – in so far as nothing to the contrary is agreed upon – ex works exclusive of freight, customs duty, import charges, packaging and VAT.

Payment

Standard payment is payment in advance. In the case of long-term business relations the following may be agreed upon:

- net for payment 30 days from date of invoice, or a 2% discount for payment within 14 days from date of invoice
- available for European customers:

 SEPA Direct Debit, i.e. payment collection on the 3rd working day of the month following the month in which invoices are due, considering the 2% cash
 ...

Prerequisite for being granted a discount shall be the settlement of all previously due and undisputed invoices. Where the due date has lapsed, all open invoices shall be immediately due for payment and the usual bank interest shall be charged.

Shipment

Delivery is for the account and at the risk of the purchaser. Shipping costs will be charged to the customer. Partial shipments initiated by us are free of charge. Costs for special shipping requests (e.g. express) will be charged for. Deliveries will be handled as fast as possible. Information regarding delivery periods, however, shall remain non-binding, in so far as they have not been agreed upon in writing and explicitly designated as binding. The delivery period shall be considered to have been complied with upon notilcation of readiness for dispatch if the shipment is not possible through no fault of the supplier. Reasona-

ble partial deliveries as well as acceptable deviations from the order quantities are permissible. The supplier shall at his discretion select packing, shipping route and method of shipment. Risk shall pass to the purchaser upon leaving the supplier's premises, even if dispatched carriage free. With regard to delays of dispatch attributable to the purchaser, risk shall pass to the purchaser upon notilcation of readiness for dispatch. In the event of Acts of God and other unforeseeable, exceptional and excusable circumstances or any disruptions of operations – even if responsibility is caused by an upstream supplier – the delivery period shall be extended by a reasonable period. Payment of damages is excluded.

arranty

Complaints to the supplier must be raised immediately in writing, at the latest 10 days after receipt of the delivery. In case of hidden deficiencies the period is extended by one week after discovery of the defect. The warranty claims shall come under the statute of limitations after one year has passed since delivery, in so far as a longer compulsory legal period is not applicable. The supplier shall not be obligated to supply a warranty if and in as far as the purchaser can invoke the limitation period of the warranty claim against his customers. With regard to substantiated complaints, the supplier shall at his discretion be obligated to repair or replace the delivery free of charge. Should the supplier not fulfil these obligations within a reasonable period of time, the purchaser shall after a reasonable extended period has elapsed without effect be entitled to reduce the purchase price or withdraw from the contract. Replaced items shall be returned to the supplier freight collect upon request. Repairs will be performed at supplier's expense. For reasons of hygiene, repairs can only be carried out on unused or laundered goods.

Replacement

For replacement only un-changed units in their original packaging and in pristine condition can be returned free of charge to the supplier. The supplier reserves the right to refuse return shipments. Also in the case of goods returned on a goodwill basis the supplier reserves the right to charge a handling fee of 20 % of the net sales value. Made to measure and custom-made products can not be replaced.

Retention of title

The deliveries shall remain the property of the supplier until all the remaining claims of the supplier against the purchaser have been fulfilled. Unpaid goods may not be pledged or transferred to others. In case of resale of the goods the resulting outstanding accounts to third parties shall be – to the amount of our credit – transferred without reservation to us.

Indemnity

In all cases in which the supplier is obliged to pay damages due to contractual or statutory regulations – even in case of unlawful acts – the supplier shall be liable only in as far as he or his legal representatives or agents are guilty of intent or gross negligence. The limitation of liability does not apply to injuries to the body or those resulting in health impairment or death that are due to gross breach of duty on the part of the supplier or gross or wilful breach of duty on the part of a legal representative or a person employed by the supplier.

Data Protection

We collect, process and use your personal data (first name, family name, address) to execute and perform the contract. Where required, your data are also transmitted to affiliated companies and/or our service partners who use the data to process the order. In addition, we use the data, subject to your consent (essential), to further nurture customer relations and to send you information regarding our current offers and prices, unless you have objected to receiving promotional material. You are able to revoke at any time in the future any consent you granted at the following address: BELSANA Medizinische Erzeugnisse, Laubanger 23, 96052 Bamberg, Germany or per email at infoexport@belsana.de

Place of performance and place of jurisdiction

As far as legally permissible, the place of performance for all mutual obligations shall be Bamberg for both parties. The place of jurisdiction shall be Bamberg, also for any proceedings in connection with documents and cheques. This agreement shall be subject to German law, with the exclusion of the UN Acquisitions Law (CISG).

Date August 2017

47

BELSANA

Should you have any questions regarding our products and services, we will be more than willing to assist you. We look forward to receiving your request. Our export team can be contacted as follows:

Phone + 49 951 30 93 03-761 Fax + 49 951 30 93 03-465 E-Mail infoexport@belsana.de

