

Genu Carezza

2




Features Benefits

- Lateral and medial coil spring
-  **Low stabilisation of the knee**
- Innovative TriTech material
-  **Very comfortable, secure fit**
- Two elastic straps for optional use
-  **Optimum fit for any leg shape**
- Anatomical shape with soft finish
-  **Good fit**

Art. no. 8358

Order no.: Article number=size | Example for ordering: 8358=M-7



Fits both left and right

Size	Lower leg circumference A		Thigh circumference B	
	(cm)	(inch)	(cm)	(inch)
XXS	29–32	11.4–12.6	36–40	14.2–15.8
XS	32–35	12.6–13.8	40–44	15.8–17.3
S	35–38	13.8–15.0	44–48	17.3–18.9
M	38–41	15.0–16.1	48–52	18.9–20.5
L	41–44	16.1–17.3	52–56	20.5–22.0
XL	44–48	17.3–18.9	56–61	22.0–24.0
XXL	48–51	18.9–20.1	61–67	24.0–26.4

Indications

- Knee pain with no or only minimal ligament laxity and with no or only slight feeling of instability
- Irritation and swelling of the knee joint
 - Tendomyopathy
 - Insertional tendinosis
 - Knee pain caused by strain
 - Soft tissue irritation following injuries and/or surgeries
 - Slight knee osteoarthritis

Effects

- Relieves and stabilises the knee joint
- Improves proprioception
- Supports sensorimotor function
- Promotes the resorption of oedemas and haematomas
- Can relieve pain

Genu Carezza, wraparound




Features Benefits

- Open design
-  **Easy to put on and take off**
- Lateral and medial coil spring
-  **Low stabilisation of the knee**
- Innovative TriTech material
-  **Very comfortable, secure fit**
- Two elastic straps for optional use
-  **Optimum fit for any leg shape**
- Anatomical shape with soft finish
-  **Good fit**

Art. no. 8362

Order no.: Article number=size | Example for ordering: 8362=L-7



Fits both left and right

Size	Lower leg circumference A		Thigh circumference B	
	(cm)	(inch)	(cm)	(inch)
XXS	29–32	11.4–12.6	36–40	14.2–15.8
XS	32–35	12.6–13.8	40–44	15.8–17.3
S	35–38	13.8–15.0	44–48	17.3–18.9
M	38–41	15.0–16.1	48–52	18.9–20.5
L	41–44	16.1–17.3	52–56	20.5–22.0
XL	44–48	17.3–18.9	56–61	22.0–24.0
XXL	48–51	18.9–20.1	61–67	24.0–26.4

Indications

- Knee pain with no or only minimal ligament laxity and with no or only slight feeling of instability
- Irritation and swelling of the knee joint
 - Tendomyopathy
 - Insertional tendinosis
 - Knee pain caused by strain
 - Soft tissue irritation following injuries and/or surgeries
 - Slight knee osteoarthritis

Effects

- Relieves and stabilises the knee joint
- Improves proprioception
- Supports sensorimotor function
- Promotes the resorption of oedemas and haematomas
- Can relieve pain