

# Genu Sensa



## Features & Benefits

Three-dimensional flat-knitted fabric with targeted compression zones

### • Optimum anatomical fit and excellent wearer comfort

Skinguard

### • Antibacterial, odour-inhibiting and skin-friendly

Patented open mesh structure in the hollow of the knee

### • Particularly soft and elastic

Lateral and medial coil spring

### • Low stabilisation of the knee

Integrated, anatomically shaped Technogel pad


### • Targeted micro-massage

## Art. no. 50K15

Order no.: Article number=size-colour | Example for ordering: 50K15=S-1

(V-Max version: Order no.: Article number=V-size-colour | Example for ordering: 50K15=V-S-1)

Fits both left and right 1=Grey 0=Beige



Size	Lower leg circumference A		Thigh circumference B		V-Max Thigh circumference B	
	(cm)	(inch)	(cm)	(inch)	(cm)	(inch)
XXS	25–28	9.8–11.0	35–38	13.8–15.0	–	–
XS	28–31	11.0–12.2	38–41	15.0–16.1	–	–
S	31–34	12.2–13.4	41–44	16.1–17.3	–	–
M	34–37	13.4–14.6	44–47	17.3–18.5	47–50	18.5–19.7
L	37–40	14.6–15.8	47–50	18.5–19.7	50–53	19.7–20.9
XL	40–43	15.8–16.9	50–53	19.7–20.9	53–56	20.9–22.0
XXL	43–46	16.9–18.1	53–56	20.9–22.0	56–59	22.0–23.2

## Indications

Knee pain with no or only minimal ligament laxity and with no or only slight feeling of instability

- Irritation and swelling of the knee joint
- Tendomyopathy
- Insertional tendinosis
- Knee pain caused by strain
- Soft tissue irritation following injuries and/or surgeries
- Knee osteoarthritis

## Effects

- Continuously optimises soft tissue compression
- Activates muscular stabilisation
- Improves proprioception
- Supports sensorimotor function
- Promotes the resorption of oedemas and haematomas
- Relieves and guides the patella
- Can relieve pain